

Performance Excellence and Easy Operation

CONTROLLED BY COMPUTER
STERILIZER
XXIII
SERIES
UDONO STERILIZATION SYSTEM



Ethylene Oxide Gas Sterilizer


udono limited

Focusing on the Fundamentals of a Sterilization System

Offering Sterilization, Purification and Safety

As a company specializing in the manufacture of sterilization equipment, Udono, while always breaking new ground, has continued supplying superior sterilization equipment required by the times.

Based on our abundant technological experience accumulated over many years, we have taken a fresh look at the fundamentals of ethylene oxide gas sterilization systems. Instead of giving the present system an excessive number of functions, we have focused on the basics, creating a system with superior sterilization capabilities and safety features.

The Color Liquid Crystal Display(LCD) Exemplifies its High Efficiency

The sterilizer, X-III was developed to provide everything: enhanced state-of-the-art technology, complete and absolute sterilization, easy to use operator-friendly equipment, preventive maintenance internal systems, an internal operating system to prevent malfunctions, and its own safety mechanism. With this piece of equipment an enlarged vision has been reach and is easily visualized.

Preconditioning Mode

For ethylene oxide gas sterilization, temperature and humidity options are the most important factors. By employing the process of preconditioning, the most desired environment for sterilization is achieved, and the failure of the attainment of sterilization, because of a temperature of humidity deficiency, is thwarted. The preheat of the sterilizing object will begin. Not only the temperature but the humidity can be adjusted. Even with just a little preconditioning time, it will greatly affect the sterilization process.

The Process Trend

Since only the temperature is not dependable, the pressure trend neutralizes this problem and instills total dependability. The large color LCD, clearly allows for the visualization of both trends - pressure and temperature, which fluctuates constantly for a maximum of 18 hours.

Aeration System



The pulse-aeration will continuously exhaust the EO gas and simultaneously intake fresh air until the EO gas has dissipated. Even after the sterilization process has been completed, if the door is not opened, and extended aeration control will be activated.

Ending Information/Reservation Timer

After the sterilization begins, the computer computes the time remaining to completion of the process. This is called Finish Information. This convenient function also reports the completion of the process. If you set the Reservation Timer to begin the sterilization process for up to 100 hours later, it will do so.

The Data Printout Recorder

A printout of the dates, starting times and completion times is absolutely essential in any sterilization process. Not only is the aforementioned recorded by the X-III but the process employed is also listed each time. If an error occurs during the printout, it also will be duly recorded.



Ethylene Oxide Gas Sterilizer

XIII SERIES

Compact Design

Although the system's capacity has been enhanced, we have succeeded in keeping the size of space required for installation at a minimum. This is achieved by giving consideration to external factors such as the space required for the piping. The installed dimensions and the size of space required for servicing are also small.

Validation

The implements of calibration and validation of the gauges attached to the sterilizer is a very effective way of ensuring the guarantee of complete and total pathogen-free sterilization. The nozzle, which is needed to measure the support temperature distribution inside the chamber, is to become standard equipment on all future Udono's sterilizers.

Error Display

The color LCD visually annunciates an error message whenever trouble occurs. Out of consideration for safety, this simple and easy to understand message allows for the problem to be directly dealt with.

The Service Display

The schedule for the drain removal of the air compressor, the exchange of hot water, and the exchange of the gas cylinder are reported. The functions are simply done and easily understood.

Excellent Safety Features

The system is equipped with a large number of safety features including the start circuit interlock (for detecting insufficient locking of doors), faulty gas supply alarm, alarms indicating pressure drop and overpressure during sterilization, safety mechanisms for preventing the mishandling of doors during operation and for preventing the two doors from being opened simultaneously, or when there is residual pressure in the chamber. Additional safety is provided by safety circuits such as those for ensuring proper sterilization temperature and for resuming operations after a power stoppage at the point at which the process was before the stoppage.

PREPARATION (Jacket heating phase)

When power is turned on water supply starts. After the proper water level is reached the heater and the water circulation pump are activated, keeping the jacket constantly warm.

VACUUM PHASE

Using a high performance vacuum pump, vacuum is attained inside the chamber to create the proper sterilization conditions. When the desired vacuum state is reached it is maintained by the system.

HUMIDITY SUPPLY AND ADJUSTMENT

To obtain the proper level of humidity necessary for sterilization, a specified amount of fresh water is supplied. Water supply is continued for a specified amount of time to ensure sufficient distribution throughout the chamber.

GAS SUPPLY AND STERILIZATION PHASE


After a specified amount of sterilization gas is supplied under pressure, the sterilization timer is activated and the sterilization process starts and continues for the specified

PURIFICATION PHASE (Vacuum/Air supply)

The system uses the pulse aeration method. A cycle consisting of vacuum attainment → vacuum maintenance → air admission → air maintenance is repeated for a specified number of times. The process uses purified air.

PROCESS COMPLETION

When all sterilization phases are completed it is indicated by the system. The safety mechanism guarding the system even after the sterilization process has been completed is an invisible feature ensuring the greater safety of the user.





udono limited

2-1-9 Motoyokoyama-cho, Hachioji-shi, Tokyo 192-0063, Japan
 Telephone +81-42-642-6301 E-mail sales@udono.com
 Facsimile +81-42-644-8350 http://www.udono.com

XIII-R

Swing Type Doors

XIII-U

Slide Type Doors



GX III - R10




GX III - R14



GX III - U55



GX III - U67



GX III - U14

STANDARD SPECIFICATIONS

Description

Ethylene Oxide Gas Sterilizer

Material:

Chamber/Door: SUS-304L

Jacket: SS-400

Performance

Sterilizing Pressure:

Maximum pressure - 0.18MPa

Sterilizing pressure - 0.12MPa

Door Gasket: Moving Gasket

- *1. #1 Vertical Sliding Door
 #2 Horizontal Sliding Door
 #3 Swing Door
- *2. Number loadable for 27 cm round type cast
- *3. S = Shelves
 MT = Loading Trolley
 SF = Semi-floor Loading Cart
 F = Floor Loading Cart

INSTALLATION REQUIREMENTS

The following are required for installation.

- *4. Steam : Stop with a valve
 Pressure supply 0.3~0.65MPa
- *5. Water : Stop with a valve
 Pressure supply 0.1~0.25MPa
- *6. Drain : Stop with a plug
 Lay-out pipe independently
- *7. Gas Exhaust (opt.): Stop with a plug
- *8. Air Compressor : Stop with a valve
 Pressure Supply 0.65 ~ 0.75MPa
- *9. Gas : Stop with high pressure valve
 Capable of pressure more than 7.0MPa
- *10. Power Source : AC200V 3 φ
 (50/60Hz)
 (AC adjustable)
 Air compressor : AC200V 3 φ
 (50/60Hz) 5A

- Use an exclusive pipe for venting out water, exhaust, and drain. Temperature may be higher than 100°C.

- Install drain pipes, steam exhaust pipes and safety valve exhaust pipes as required.

- For floor loading sterilizer, space is required below surface for installation (floor level - 150 mm).

System specifications		Chamber dimensions			Capacity	Door type	Cast	Loading	System outer dimensions										
Type		W	H	D(mm)	(ℓ)	★1	★2	System ★3	W	H	D								
GXⅢ-U556-S		500 × 500 × 600			150	Slide#1	4	S/MT	1,150	1,800	1,100								
559-S		500 × 500 × 900			230		6		1,150		1,400								
677-S		670 × 670 × 700			310		8		1,220		1,200								
6710-S		670 × 670 ×1000			440		12		1,220		1,500								
6713-S		670 × 670 ×1300			580		16		1,220		1,800								
1410-S		670 ×1400 ×1000			940	Slide#2	24	SF/F	2,200		1,500								
1412-S		670 ×1400 ×1200			1,120		32		2,200		1,700								
559-D		500 × 500 × 900			230	Slide#1	6	S/MT	1,150	1,800	1,210								
6710-D		670 × 670 ×1000			440		12		1,220		1,310								
6713-D		670 × 670 ×1300			580		16		1,220		1,610								
1410-D		670 ×1400 ×1000			940	Slide#2	24	SF/F	2,200		1,450								
1412-D		670 ×1400 ×1200			1,120		32		2,200		1,650								
GXⅢ-R107-S		670 ×1000 × 700			470	Swing#3	12	S/MT	1,370	1,800	1,100								
1010-S		670 ×1000 ×1000			670		18		1,370		1,400								
1013-S		670 ×1000 ×1300			870		24		1,370		1,700								
1410-S		670 ×1400 ×1000			940		24	SF/F	1,450		1,500								
1412-S		670 ×1400 ×1200			1,120		32		1,450		1,700								
1010-D		670 ×1000 ×1000			670		18	S/MT	1,370		1,260								
1013-D		670 ×1000 ×1300			870		24		1,370		1,560								
1410-D		670 ×1400 ×1000			940		24	SF/F	1,450		1,290								
1412-D		670 ×1400 ×1200			1,120		32		1,450		1,490								
Installation requirements		Steam ★4		Water supply ★5		Drain ★6	Gas Exhaust ★7	Air pressure ★8	Gas ★9	Power supply ★10									
Type		kg/hr	Caliber	ℓ/min	Caliber	Caliber	Caliber	Caliber	Caliber	Ampare									
GXⅢ-U556-S		5	20A	8	20A	40A	25A	15A	8mm	30A	15A								
559-S											10	20A							
677-S																			
6710-S		10									20A								
6713-S																			
1410-S		25		10		50A				30A	30A								
1412-S																			
559-D		5		8		40A				20A	20A								
6710-D		10																	
6713-D		15																	
1410-D		25		10		50A				30A	30A								
1412-D																			
GXⅢ-R107-S		15	20A	8	20A	40A	25A	15A	8mm	25A	25A								
1010-S																			
1013-S																			
1410-S		20				50A				30A	30A								
1412-S																			
1010-D		15		8		40A				25A	25A								
1013-D		20																	
1410-D		25		10		50A				30A	30A								
1412-D																			

S : Single door D : Double door

Subject to modification without notice to improve quality and performance.